

Data Management for Mine Labs

ASSAYNET

Release Notes

Sample Track 1.2.4



© Copyright 2024 Datamine Software

All Rights Reserved Confidential and Proprietary

Published: Thursday, March 14, 2024

Legal Disclaimer

The product described in this documentation may be connected to, and/or communicate information and data via, a network interface, which should be connected to a secure network. It is your sole responsibility to ensure a secure connection to the network and to establish and maintain appropriate measures (such as but not limited to the installation of firewalls, application of authentication measures, encryption of data, installation of antivirus programs, etc.) to protect the product, the network, your systems, and the interface against any kind of security breach, unauthorised access, interference, intrusion, leakage, damage, or corruption or theft of data. We are not liable for damages or losses related to any such security breach, unauthorised access, interference, intrusion, leakage, damage, or corruption or theft of data.



Contents

Overview	4
Sample Track Release Notes	5
Version 1.2.4 (8-Mar-24)	5
Version 1.2.3 (29-Feb-24)	5
Version 1.2.2 (15-Feb-24)	5
Version 1.2.1 (13-Feb-24)	5
Version 1.2.0 (19-Jan-24)	5
Version 1.1.1 (17-Nov-23)	5
Version 1.1.0 (19-Jun-23)	6
Version 1.0.3 (14-Dec-22)	6
Version 1.0.2 (01-Dec-22)	6



Overview

AssayNet is a client/server laboratory information management system (LIMS) built specifically for environmental and minerals laboratories. Robust and highly scalable, the program was written to conform to the American Society for Testing Materials (ASTM) level II specification for LIMS. AssayNet contains all the features you would expect in an industrial-strength LIMS. From sample reception, logging and preparation through analysis and reporting, the program has been designed to work the way you do. Its intuitive, menu-driven interface is easy to learn, yet feature rich, providing access to a complete quality control (QC) program, quotations and invoicing, stock levels and inventories, sample storage and more.

Further Information

This document includes cumulative release notes for AssayNet Sample Track.

Release notes for other versions of AssayNet are available via the Support Portal <https://www.dataminesoftware.com/support/>.



Sample Track Release Notes

Version 1.2.4 (8-Mar-24)

- Add submittal name filter.
- Add confirmation message when manually entering a new sample ID.
- Update app color scheme.
- Automatically connect when scanning auto-login barcode.

Version 1.2.3 (29-Feb-24)

- Add option to open sample details after scanning.
- Fix submittal sample count.

Version 1.2.2 (15-Feb-24)

- Fix scanning barcodes with surrounding asterisks.

Version 1.2.1 (13-Feb-24)

- Add show password option.
- Allow deleting selected submittals.
- Fix app not working with newer Zebra software version.

Version 1.2.0 (19-Jan-24)

- Add mapping feature to show planned and captured sample GPS locations.

Version 1.1.1 (17-Nov-23)

- Check DB support for GPSCoordinates before calling procedures.



Version 1.1.0 (19-Jun-23)

- Add support for scanning series samples into submittals.
- Fix sample numbers not updating in submittal list.

Version 1.0.3 (14-Dec-22)

- Add better logging and more detail to log export.

Version 1.0.2 (01-Dec-22)

- Support scanning barcodes with either Sample ID strings or Submittal Sample ID e.g *SS000001234*.



Datamine enables efficient and sustainable mining through the application of world-leading technology and services.

Read the Docs

docs.dataminesoftware.com

Get in Touch

sales@dataminesoftware.com

www.dataminesoftware.com/support



Find Us

AUSTRALIA | BRAZIL | CANADA | CHILE | CHINA |
ECUADOR | GHANA | INDIA | INDONESIA | KAZAKHSTAN |
MALAYSIA | MEXICO | MONGOLIA | PERU | PHILIPPINES |
SOUTH AFRICA | TURKEY | UNITED KINGDOM | USA |
UZBEKISTAN



www.dataminesoftware.com



ASSAYNET



DATAMINE